

DIGITAL IMAGING PLATES >>>

Imaging plates for computed radiography application

Phosphor plates can be erased and reused thousands of times and will provide faster image access. Utilization of computed radiography also allows for manipulation of the image either in the lab or field, and further enables you to electronically share these images.

Product Features

- High absorption efficiency
- New electron beam cured topcoat creates high tolerance to mechanical wear
- Smooth surface for improved signal-to-noise ratio
- All sizes and customization available
- Flexible, portable, or rigid cassette application

Systems Supported

- Agfa
- Carestream (Kodak)
- Fuji
- Orex
- Radlink
- Allpro
- Lumisys/ACR 2000
- Alara
- Konica

Industries Supported

- Medical
- Veterinary
- Chiropractic
- Dental
- NDT
- Podiatry

Cleaning/Handling of CR Plates

- It's extremely important to use the correct cleaning solutions on imaging plates. Using the wrong one can permanently damage the plate and/or compromise image quality.
- Wear gloves when handling CR cassettes and plates.
- Contact Rochester Cassette for the correct cleaning solution for your CR cassettes and imaging plates.

WIDE VARIETY STOCKED FOR IMMEDIATE SHIPMENT



All imaging plates are guaranteed to be free of manufacturer defects

ISO 9001:2000 CERTIFIED

No surprises – all work is guaranteed



Rochester Cassette
Sales & Service

50 Old Hojack Lane, Hilton, New York 14468

1.800.528.6311 • 585.392.4336

585.392.5575 (fax)

www.rochestercassette.com

sales@rochestercassette.com